

## Total viable counts in different foods: Test tube method versus Inlabtec for serial dilutions

Sample code	Sample type	logcfu/g standard	log cfu/g inlabtec
1 R	yogurt drink	7.87	7.89
3 R	fish pate	5.66	5.89
4 R	vegetable juice	2.78	2.60
5 R	salad with protein	6.17	6.30
6 R	fresh parsley	6.41	6.51
7 R	pork mince	3.74	3.71
8 R	chicken fillets	2.00	2.48
10 R	fresh prawns	3.20	3.26
3 F	fish pate	5.79	5.83
4 F	vegetable juice	2.85	3.08
5 F	salad with protein	6.30	6.29
6 F	fresh parsley	6.47	6.32
7 F	pork mince	3.91	3.95
8 F	chicken fillets	2.48	2.48
22	Flour	4.95	5.04
25	Flour	4.40	4.65
27	Flour	5.06	5.18
28	Flour	4.23	4.43
29	Flour	4.35	4.67
30	Flour	4.71	4.72
31	Flour	5.32	5.43
32	Flour	4.00	4.18
33	Flour	4.80	5.04
34	Flour	3.56	3.76
18	chicken noodle pot	9.30	9.26
19	Egg and tomato salad	3.24	3.45
64	Grated Cheese	7.04	7.34
67	Blue Cheese	8.04	8.33
236	Jam Spong	5.15	5.53
327	Beef lasagne	8.76	8.89
344	Tomato soup	5.71	5.89
350	Cheese and tomato salad	8.90	9.00
354	Chicken noodle salad	8.48	8.62
357	Mozerella and tomato salad	9.26	9.34

Correlation Coefficient: 0.998

